**CERTEST FOB CARD**  
One step Fecal Occult Blood Card Test

**INTENDED USE**

The CerTest FOB Card is a one step coloured chromatographic immunoassay for the qualitative detection of human hemoglobin in stool samples.

**INTRODUCTION**

Colorectal cancer is a leading cause of illness and death in Western world. Screening with fecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps, will bleed and may be detected by an occult blood test. When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact hemoglobin, intact heme, and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut. Immunochemical tests detect intact or nearly intact human hemoglobin, being a very specific technique for detecting loss of blood from the lower intestine, because blood from lower sites is less degraded during transit.

**PRINCIPLE OF THE TEST**

The CerTest FOB Card is a qualitative immunochromatographic assay for the determination of human hemoglobin in stool samples. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, against human hemoglobin. During testing, the sample is allowed to react with the coloured conjugate (anti-human hemoglobin mouse monoclonal antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

**STORAGE AND STABILITY**

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

**PRECAUTIONS**


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

**SPECIMEN COLLECTION AND PREPARATION**

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

**Specimen preparation (see illustration):**

- (1) Unscrew the tap and use the stick to pick up a little sample, if the stool sample was liquid take 100 µL using a pipette, and add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

**Precautions:** Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine, or if they have strained during bowel movement.

**MATERIALS PROVIDED**

- Device tests
- Instructions for use
- Stool collection tubes-sample diluent


**MATERIALS REQUIRED BUT NO PROVIDED**

- Specimen collection container
- Disposable gloves
- Timer


**TEST PROCEDURE**

**Allow the test, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.**

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the top (3).
2. Remove the CerTest FOB Card device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense exactly 5 drops or 150 µL into the circular window marked with an arrow (4).
- 4.- Read the result at **5 minutes** (the coloured bands appear). Do not read a test result after more than 10 minutes.

**INTERPRETATION OF RESULTS (PLEASE REFER TO THE ILLUSTRATION BELOW)**

**NEGATIVE:** Only one GREEN band appears across the central window in the site marked with the letter C (control line).

**POSITIVE:** In addition to the GREEN control band, a RED band (test line) also appears in the site marked with the letter T (result line).

**INVALID:** a total absence of the control coloured band (GREEN) regardless the appearance or not of the result line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the tests with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

**NOTES ON THE INTERPRETATION OF RESULTS**

The intensity of the RED coloured band in the result line region (T) will vary depending on the concentration of human hemoglobin in the specimen. However, neither the quantitative value, nor the rate of increase in hemoglobin can be determined by this qualitative test.

**QUALITY CONTROL**

Internal procedural controls are included in the test. A GREEN line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

**LIMITATIONS**

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear or absence of the control coloured band).
3. The intensity of the lines may vary from very strong at high hemoglobin concentrations to faint when the hemoglobin concentration is close to the sensitivity limit of the test.
4. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine, or if they have strained during bowel movement.
5. Positive results confirm the presence of human hemoglobin in fecal samples; nevertheless, it can be due to several causes, besides colorectal bleeding, such as haemorrhoid, blood in urine or stomach irritations. A positive results should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood in the stool.
6. Negative results do not exclude bleeding since some polyps and colorectal cancers may bleed intermittently or not at all during certain stages of the disease. Additionally, blood may not be uniformly distributed in stool samples.
7. This test may be less sensitive for detecting upper gastrointestinal bleeding because blood degrades as it passes through the gastrointestinal track.

**PERFORMANCE**

**Sensitivity**

A sample containing human hemoglobin at concentration equal to or higher than 50 ng/mL produces a positive result using CerTest FOB Test. In some cases sample containing human haemoglobin at concentrations less than 50 ng/ml can also be tested as positive.

Different hemoglobin dilutions were tested directly in the extraction buffer or spiked in a negative stool sample in accordance with the kit instructions.

The detection of human hemoglobin with CerTest FOB Test showed >99% of sensitivity compared with a commercial guaiac assay.

**Specificity**

The CerTest FOB Test is specific for human hemoglobin and do not show any cross reaction with hemoglobin from bovine and pig.

The detection of human haemoglobin with CerTest FOB Test showed >99% of specificity compared with a commercial guaiac assay. The use of mouse monoclonal antibodies in the elaboration of CerTest FOB Test assures high degree of specificity for the detection of human hemoglobin.

**REFERENCES**

1. Towler BP, Iwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

**SYMBOLS FOR IVD COMPONENTS AND REAGENTS**

	In vitro diagnostic device		Lot code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Do not use if package damaged

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1. Towler BP, Inwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Faecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

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